

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to www.agilent.com/chem.



Structural Elucidation of Ginsenosides Using the Varian 500-MS Ion Trap Mass Spectrometer

Mahalakshmi Rudrabhatla
Varian, Inc.

Introduction

Ginsenosides are active components of *Panax ginseng*¹. Among the 30 ginsenosides identified so far, ginsenosides Rb1, Rb2, Rc, Re and Rd are the most abundant and are associated with the pharmacological activity of the plant material². Since ancient times, traditional Chinese medicine (TCM) has been used in the treatment of atherosclerosis, arthritis, asthma, diabetes, stroke, multiple sclerosis and endotoxin liver injury³. Ginseng root (*Panax ginseng*) is the major ingredient in TCM. Active components of ginseng root are ginsenosides that fall under the category of nutraceuticals (food components/products that have therapeutic and pharmacological properties). With an estimated 15 million patients at risk of potentially adverse drug-herb interactions, there is a renewed interest in the analysis of these compounds.

The main problem in the analysis of ginsenosides is their structural similarity. They are described as glycosides consisting of an aglycone moiety that is typically a triterpenoid or steroid and one or more covalently linked sugar monomers. Most ginsenosides contain multiple oligosaccharide chains at different positions. Analysis of these compounds can therefore be quite complicated. ESI-TOF, triple quadrupole MS/MS and Ion trap MSⁿ methods have been developed so far for the structural elucidation and analysis of these compounds in complex matrices. Among these methods, MSⁿ analysis in an ion trap mass spectrometer permits multiple isolation and fragmentation stages, ensuring that product ions in each stage are specifically related to the precursor ion from that particular stage. It is possible to confirm the authenticity of the pharmaceutical ginseng products and differentiate between their active ingredients using ion trap MSⁿ fragmentation patterns.

This application note demonstrates the use of an LC/MSⁿ method for the structural elucidation of ginsenosides Rb1, Rb2, Rc, Re and Rd on the Varian 500-MS Ion Trap LC/MS. This method was validated in Korean ginseng extract. The information obtained from the ion trap MSⁿ and MRM could be used for determining the structure of ingredients and for quality control purposes in complex matrices.

Instrumentation

- Varian 500-MS LC Ion Trap with an ESI source
- Varian 212-LC Binary Gradient LC/MS Pumps
- HTS Pal AutoSampler

Materials and Reagents

Ginsenosides Re, Rc, Rd, Rb1 and Rb2 were obtained from Chromadex Inc, Santa Ana, CA. HPLC grade water and acetonitrile were obtained from Fisher Scientific.

Sample Preparation

Powdered Korean ginseng (0.3 g) was ultra sonicated in 6 mL methanol for 30 minutes. The extract was filtered initially using Whatmann #1 filter paper. This was further filtered using a 0.2 µm syringe filter.

LC Conditions

Column: Polaris™ C18-A 3 µm, 150 x 2 mm ID
(Varian Part Number A2001150X020)

Buffer A: Water

Buffer B: Acetonitrile

Injection Volume: 50 µL

LC Program:	Time (min:sec)	%A	%B	Flow (µL/min)
	0:00	75.0	25.0	200.0
	40:00	60.0	40.0	200.0
	40:01	75.0	25.0	200.0
	50:00	75.0	25.0	200.0

API Conditions

Ionization Mode: ESI Positive

API Drying Gas: 15 psi at 350 °C

API Nebulizing Gas: 25 psi

Needle: 4000 V

Shield: 600 V

Detector: 1400 V

TurboDDS™ Scan Conditions

Scan Type: TurboDDS

Survey Scan Start Mass: 900 *m/z*

Survey Scan End Mass: 1200 *m/z*

RF Loading: 100%

MSⁿ Depth: *n* = 3

Auto Mass Range: Yes

Isolation Window: 3.0 *m/z*

An exclusive include list containing molecular masses of all five ginsenosides and their respective MS² and MS³ fragments was created. TurboDDS™ was ordered to trigger the masses given in the include list based on abundance. The exclusive include abundance list for all five ginsenosides is given in Table 1.

Table 1. Exclusive include abundance list for all five ginsenosides.

Ginsenosides	MS ¹	MS ²	MS ³
Re	967.0– 972.0	787.0 – 791.0	347.0 – 350.0
Rb1	1130.0 – 1133.0	787.0 – 791.0 363.0 – 366.0	363.0 – 366.0
Rc/Rb2	1100.0 – 1102.0	787.0 – 790.0 333.0 – 336.0	363.0 – 366.0
Rd	967.0 – 972.0	787.0 – 791.0	363.0 – 366.0

Results & Discussion

Initially, all five ginsenosides (Rc, Re, Rb1, Rb2 and Rd) were infused at 300 pg/μL concentrations on the Varian 500-MS. MS/MS and MS³ data was acquired automatically using the TurboDDS™ software in a data dependent fashion to investigate in detail the structural elucidation of these ginsenosides. Even though ginsenosides ionize well in the negative ion mode under normal electrospray conditions, structural elucidation cannot be done at the MS² level since all five ginsenosides are similar. In the positive ESI mode protonated sodiated adduct ions [M+Na]⁺ are formed when high drying gas temperatures are applied (350 °C). MS² and MS³ analysis of these ions gave much clearer fragmentation, thus facilitating structural elucidation.

The fragmentation pattern of ginsenoside Re on the 500-MS using the TurboDDS software is given in Figure 1. The main ion observed for ginsenoside Re in the positive mode is [M+H] + Na = 969.2. MS² fragmentation yielded a single product ion at *m/z* = 789.1 that was produced by the loss of a glucose molecule from the sodiated molecular ion 969.2. At the MS³ level, the MS² product ion 789.1 is cleaved into two parts. The detected part at MS³ level is *m/z* = 348.8, which comes from the cleaved disaccharide moiety.





Figure 1. Ginsenoside Re fragmentation on 500-MS (left). Structure of ginsenoside Re (right).





Figure 2. Ginsenoside Rd fragmentation pattern on 500-MS (left). Structure of ginsenoside Rd (right).

Ginsenoside Rd has the same sodiated protonated molecular ion (*m/z* = 969.1) and MS² product ion (*m/z* = 789.2) as ginsenoside Re. Interestingly, ginsenoside Rd could be differentiated from ginsenoside Re by the fragmentation at MS³ stage on 500-MS. Fragmentation of ginsenoside Rd yielded a disaccharide ion at *m/z* 364.9 at the MS³ stage (Figures 1 and 2).

The main precursor ion observed for ginsenoside Rb1 is at *m/z* 1131.3, and is subsequently fragmented in the MS/MS experiment into two product ions at *m/z* 789.3 and *m/z* 364.9 (Figure 3). At the MS³ level, *m/z* 789.3 is further fragmented to yield a disaccharide ion at *m/z* 364.8.

Isomeric ginsenosides Rc and Rb2 have the same fragmentation patterns. Two fragments were obtained at the MS² level, one at *m/z* = 789.2 and the other at 334.9 for the arabinose disaccharide. The other disaccharide *m/z* 364.8 is cleaved off from the ion 789.2 at the MS³ level (Figure 4).

All five ginsenosides analyzed have the same MS² level main product ion corresponding to *m/z* = 789.1. This provides a structural template that supports the fragmentation patterns proposed for these compounds.

As the saccharide molecules connected to the triterpene structures of ginsenosides are cleaved off in a characteristic pattern with the ion trap MS/MS and MSⁿ, it is possible to differentiate ginsenosides in ginseng extract. If the specific fragmentation pattern of ginsenosides is known, it helps to detect the presence of ginsenosides and specifically ginseng root in complex natural product extracts.




Figure 3. Ginsenoside Rb1 fragmentation pattern on 500-MS (left). Structure of ginsenoside Rb1 (right).




Figure 4. Ginsenoside Rc/Rb2 fragmentation pattern on 500-MS (left). Structure of ginsenoside Rc/Rb2 (right).

In the next step, ginsenosides Re, Rc, Rb1, Rb2 and Rd were separated chromatographically using a Varian Polaris™ C18-A 150 x 2.0 mm ID column. Product ions obtained by the fragmentation of sodium adducts in the positive ion mode were used for monitoring MRM transitions. Excellent separation with good resolution was obtained for all five ginsenosides (Figure 5).

Chromatographic retention times of ginsenosides and their respective MRM transitions are given in Table 2. As seen from the table, chromatographic retention times along with MS² and MS³ product ions specific for each ginsenoside could be used for generating qualitative fingerprints of these compounds in complex extracts.

Table 2. MRM transitions and LC retention times of ginsenosides studied.

Name	Retention Time (min)	Parent Ion	MS ² Product Ions	MS ³ Product Ions
Ginsenoside Re	4.751	969.2	789.1	348.8
Ginsenoside Rb1	18.485	1131.3	789.3, 364.9	364.8
Ginsenoside Rc	20.082	1101.2	789.2, 334.9	364.8
Ginsenoside Rb2	22.086	1101.2	789.2, 334.9	364.8
Ginsenoside Rd	26.317	969.1	789.2	364.9




Figure 5. Extracted ion MRM chromatograms of five ginsenosides.

In the next step, Korean ginseng extract was tested for the presence of these ginsenosides. 0.3 g of Korean ginseng powder was extracted into 6 mL methanol. This extract was diluted 500 times with 25% acetonitrile before injection on the 500-MS. As shown in Figure 6, all five ginsenosides were detected in the extract. MS² and MS³ spectra extracted from the peaks further confirmed their presence thus generating their fingerprint.




Figure 6. Five ginsenosides observed in the Korean ginseng extract.

Conclusion

This application note successfully demonstrates the use of the Varian 500-MS Ion Trap mass spectrometer with its powerful TurboDDS™ software as an ideal tool for the structural elucidation of natural products in complex natural extracts. Complex structures of five ginsenosides, which are the main components of natural extracts, could be elucidated by the interpretation of the MS/MS and MSⁿ ion trap data obtained. Detailed knowledge of fragmentation data on an ion trap could be used for the quality control of complex natural extracts or pharmaceutical products by ion trap MSⁿ when applied to structurally similar active ingredients in natural products.

References

1. Attele, A.; Wu, J.A.; Yuan, C.S. *Biochem Pharmacol.* 1999, 58:1685-1693.
2. Shibata, S.; Tanaka, O.; Shoji, J.; Saito, H. *Economic and Medicinal Plant Research.* 1985, Vol 1.
3. Kim, Y.C.; Kim, S.R.; Markelonis, G.J.; Oh, T.H. *J Neurosci. Res.* 1998, 53, 1998, 426-32.
4. Wang, X.; Sakuma, T.; Asafu-Adjaye, E.; Shiu, G.K. *Anal. Chem.* 1999, 71, 1579-1584.

Polaris, TurboDDS, Varian and the Varian logo are registered or unregistered trademarks of Varian, Inc. in the U.S. and other countries.
© 2007-2008 Varian, Inc.

These data represent typical results.
For further information, contact your local Varian Sales Office.

Varian, Inc.
www.varianinc.com
North America: 800.926.3000 – 925.939.2400
Europe *The Netherlands*: 31.118.67.1000
Asia Pacific *Australia*: 613.9560.7133
Latin America *Brazil*: 55.11.3845.0444