

Tips to Help Maximize Resolution

Mark Powell
Columns and Supplies Technical Support
August 16, 2022



What is My Method Development Plan?

1. Smaller particles and superficially porous particles offer fast, efficient analysis
2. C18 column, a general-purpose column choice
3. Simple mobile phase
 - a) Formic acid or other additive in aqueous portion (buffer salts if necessary)
 - b) Acetonitrile or methanol as organic modifier
4. Start with linear gradient (5% organic to 95% organic) for reversed-phase methods
5. Adjust mobile phase to get the desired retention and resolution
 - a) Adequate resolution of all peaks, $R_s \geq 2.0$
 - b) Retention of first peak at least $k=1$
 - c) Fastest analysis time with required resolution

Speed up method development by using shorter columns with small particle sizes. Columns like these can provide increased efficiency and resolution in a shorter time.

What Column Do I Choose?

Smaller particle size offers

- Higher efficiency, shorter column, faster method
- Increased resolution
- Better sensitivity


But make sure to consider pressure limit of instrument

Smaller diameter means


- Solvent savings

But this depends on instrument configuration and plumbing

- Bonded phase choices
 - Alternate selectivity
 - Match to pH of mobile phase
 - More robust column life




What Particle Do I Choose?



Totally Porous Particle

ZORBAX Eclipse Plus C18
4.6 x 250 mm, 5 μ m

Run time: 35 min



Poroshell Particle

InfinityLab Poroshell 120 EC-C18
4.6 x 100 mm, 2.7 μ m

Run time: 9 min

A: 0.1% Formic Acid in water, B: ACN
Gradient: 8–33% ACN in 30 or 8 min
1 or 2 mL/min, 25 $^{\circ}$ C, 254 nm
Agilent App Note, 5990-5572EN

Poroshell Particles

SPP particle	For	Maximum pressure	Typical pressure	Efficiency	Target system
1.9 μm	Highest UHPLC performance	1300 bar	Similar to sub-2 μm totally porous	~120% of sub-2 μm totally porous	1290 Infinity II
2.7 μm	UHPLC performance at lower pressures	600 bar / 1000 bar	50% of sub-2 μm totally porous	~90% of sub-2 μm totally porous	1290 Infinity II 1260 Infinity II
4 μm	Improved HPLC performance	600 bar	Typically < 200 bar	~200% of 5 μm totally porous	1260 Infinity II VL 1220 Infinity II (VL)


Column length	Recommended Use
50	High speed
100	High resolution
≥ 150	Ultrahigh resolution

Particle size	id	Optimum flow
1.9 μm	2.1 mm	0.4 – 0.5 mL/min
	3.0 mm	0.8 – 1 mL/min
2.7 μm	2.1 mm	0.4 – 0.5 mL/min
	3.0 mm	0.8 – 1 mL/min
	4.6 mm	1.5 – 2 mL/min
4 μm	3.0 mm	0.5 – 0.75 mL/min
	4.6 mm	1 – 1.25 mL/min

Factors That Affect Resolution

$$R_s = \left(\frac{1}{4}\right) N^{0.5} \left(\frac{\alpha - 1}{\alpha}\right) \left(\frac{k}{1 + k}\right)$$

Resolution Efficiency Selectivity Retention



Selectivity impacts resolution the most

- Change bonded phase
- Change mobile phase

Evaluate Different Bonded Phases

- Bonded phase affects selectivity (alpha)
- Different interactions for polar and nonpolar compounds.
- Exploit other interactions with bonded phase (e.g., pi-pi)
- Changing the bonded phase can improve selectivity/resolution,
- May reduce analysis time
- Having different bonded phases available on the same particle makes development easier





Evaluating different bonded phase chemistries early can save time in optimization and generate a more robust method

The InfinityLab Poroshell 120 Portfolio


Agilent Poroshell columns are designed for multiple separation modes

Best all around	Best for low pH mobile phases	Best for high pH mobile phases	Best for alternative selectivity	Best for more polar analytes	Chiral
EC-C18 ^A 1.9 μm, 2.7 μm, 4 μm	SB-C18 ^A 1.9 μm, 2.7 μm, 4 μm	HPH-C18 ^A 1.9 μm, 2.7 μm, 4 μm	Bonus-RP ^{A,B} 2.7 μm	SB-Aq ^{A,B} 1.9 μm, 2.7 μm, 4 μm	Chiral-V ^{A,C,D} 2.7 μm
EC-C8 ^A 1.9 μm, 2.7 μm, 4 μm	SB-C8 ^A 2.7 μm	HPH-C8 ^A 2.7 μm, 4 μm	PFP ^{A,B,D} 1.9 μm, 2.7 μm, 4 μm	EC-CN ^{A,B,C,D} 2.7 μm	Chiral-T ^{A,C,D} 2.7 μm
Phenyl-Hexyl ^A 1.9 μm, 2.7 μm, 4 μm			CS-C18 ^A ← 2.7 μm →	HILIC ^{C,D,E} 1.9 μm, 2.7 μm, 4 μm	Chiral-CD ^{A,C,D} 2.7 μm
Legend A reversed phase B can be operated at 100% aqueous C Normal phase D SFC E HILIC				HILIC-Z ^{C,D,E} 1.9 μm, 2.7 μm, 4 μm	Chiral-CF ^{A,C,D} 2.7 μm
				HILIC- OH5 ^{C,D,E} 2.7 μm	

What C18 Bonded Phase?

InfinityLab Poroshell 120	Chemistry	Pore Size	Endcapped	Carbon Load	Surface Area	Best For
EC-C18 1.9 μm , 2.7 μm , 4 μm		120 Å	Yes	10%	130 m ² /g	General Purpose Excellent peak shape and efficiency for acids, bases, neutrals
SB-C18 1.9 μm , 2.7 μm , 4 μm		120 Å	No	9%	130 m ² /g	Low pH Excellent stability and peak shape in highly acidic conditions
HPH-C18 1.9 μm , 2.7 μm , 4 μm		100 Å	Yes	Proprietary	95 m ² /g	High pH Robust performance and long lifetimes
CS-C18 2.7 μm		100 Å	Yes	Proprietary	95 m ² /g	Alternate selectivity Improved peak shape and sample capacity for basic compounds with low ionic strength mobile phases High pH

Alternative Selectivity with InfinityLab Poroshell 120 C18s




1	Ciprofloxacin	-----
2	Oxytetracycline	-----
3	Tetracycline	-----
4	Enrofloxacin	-----
5	Sulfamerazine	-----
6	Sulfamethazine	-----
7	Erythromycin	-----
8	Penicillin-G	-----
9	Oxacillin	-----

Method parameters:
 A: 0.1% formic acid in water
 B: acetonitrile
 0.4 mL/min, 0-95% B in 15 min
 0.05 µL injection
 Sample: 0.1 mg/mL in water
 Column: 30 °C, 2.1 x 100 mm, 2.7 µm
 Detection: LC/MS, ESI+, dMRM

Agilent application note: [5994-2358EN](#)

Selectivity Differences Across InfinityLab Poroshell Bonded Phases



1. Hydrocortisone 2. β -Estradiol 3. Androstatriene-3,17-dione 4. Testosterone
5. Ethinyl estradiol 6. Estrone 7. Norethindrone acetate 8. Progesterone


40–80 % Methanol in 14 min, DAD 260, 80 nm 0.4 mL/min,
2.1 x 100 mm column, 40 C, 0.1% formic acid in water and
methanol, Agilent 1260 Method Development Solution

Agilent InfinityLab Poroshell 120 CS-C18

Mobile phase pH is a method development tool for separating ionizable compounds

- With reversed-phase, ionizable analytes are more retained in their neutral state
- **Acids** are more retained at low pH
- **Bases** are more retained at high pH
- **Neutrals** are not affected by mobile phase pH

5-95% CH₃CN in 10 min, 4 min postrun, mobile phase A varies, 0.4 mL/min, 2.1 x 100 mm, 2.7 μm Agilent InfinityLab Poroshell 120 CS-C18, 30 °C, DAD: 254 nm, 80 Hz; Sample: uracil, amitriptyline, butylparaben, dipropyl phthalate, acenaphthene



Agilent application note: 5994-2274EN

What Mobile Phase Modifiers Should I Try?

Mobile Phase	Useable pH range	Recommended for Silica-Based LC Columns?	Recommended for LC/MS Use?
TFA	<1.5	Limited	No
Phosphate	1.1–3.1	Limited	No
Formic acid	<2.8	Yes	Yes
Acetic acid	<3.8	Yes	Yes
Formate	2-8–4.8	Yes	Yes
Acetate	3.8–5.8	Yes	Yes
Carbonate	5.4–7.4	Yes	Yes
Phosphate	6.2–8.2	Limited	No
Bicarbonate	6.6–8.6	Limited	Yes
Ammonia	8.2–10.2	Limited	Yes
Phosphate	11.3–13.3	Limited	No


What Organic Solvent Should I Use?

Try both

- ACN and MeOH are readily available
- Works on any bonded phase – optimize separation no matter the column choice


MeOH – Higher pressure, may give better peak shape with bases, protic solvent

Acetonitrile – Aprotic, wider UV window, stronger than MeOH



InfinityLab Quick Connect and Quick Turn Fittings

- Spring loaded design
- Easy, **no tools needed**
- Works for all column types
- Reusable
- Consistent ZDV connection



Quick Connect fitting

- Finger tight up to 1300 bar
- Hand tighten the nut, then depress the lever

Quick Turn fitting

- Finger tight up to 400 bar
- Up to 1300 bar with a wrench
- Compact design



Tips for Robust Methods

- Always start method development with a new column
- Select columns with robust properties at pH of method
- Choose a quality column with long lifetimes
- Consider batch-to-batch reproducibility
- Consider scalability of particle sizes and chemistries for downstream method transfer
- Make sure that mobile phase preparation is documented and transferrable

Agilent employs end-to-end process control for quality LC columns

www.agilent.com/chem/qualitylc



What Should I Do with a New Column?

Performance report

SERIAL NUMBER: USDAZ01333


PART NUMBER: 959758-902
COLUMN TYPE: ZORBAX RRHD Eclipse Plus C18 2.1 x 100 mm, 1.8 µm
PACKING LOT #: B09089

TEST CONDITIONS

MOBILE PHASE = 60% Acetonitrile / 40% Water
COLUMN PRESSURE = 517.2 Bar
COLUMN FLOW = 0.50 ml / min
LINEAR VELOCITY = 0.436 cm / sec
TEMPERATURE = AMBIENT (Nominally 23 °C)
INJECTION VOLUME = 1 µl

QUALITY CONTROL PERFORMANCE RESULTS FOR NAPHTHALENE

TEST VALUES	SPECIFICATIONS
THEORETICAL PLATES = 22337	MIN = 21000
SELECTIVITY = 1.90	RANGE = 1.82 - 1.92
USP TAILING FACTOR = 1.08 (@ 5% Peak Height)	RANGE = 0.98 - 1.20
k' = 4.58	



Sample components with concentrations diluted in mobile phase in the following elution order.

Peak #	Conc (ug/ml)	Sample Component
1	10	Uracil
2	400	Phenol
3	50	4-Chloro Nitrobenzene
4	80	Naphthalene


A manufacturing test chromatogram is done on a modified LC system to minimize ECV and will differ from a typical lab HPLC.

- Don't expect to get the exact same result as the performance report
- Test column performance on your instrument to have as a reference

What Should I Do With a New Column?

Benchmark new column on your system

1. Standard mix; test mix (5188-6529, 01080-68704; QC reference material)
2. Criteria like retention time, peak area, peak tailing, resolution, response, and system pressure
3. Theoretical plates
 - Monitor column over time
 - Troubleshoot



Chromatographic conditions

Sample: RRLC Checkout sample
(p/n 5188-6529)
Column: Agilent Poroshell 120
EC C18, 3 mm × 50 mm,
2.7 μm
Mobile phase: A = Water
B = Acetonitrile
Gradient: 0 min 20% B
8 min 80% B
Flow rate: 1.2 mL/min
Stop time: 8 min
Post time: 4 min
Injection volume: 1 μL
Column temperature: 30 °C
DAD: 245/10 nm
Ref 400/100 nm
Flow cell: 10 mm
Peak width: <0.025 min (10 Hz)

Mobile Phase Preparation


Specified volume ACN added to a 1 L volumetric and made to volume with H₂O

≠

Specified volume H₂O added to a 1 L volumetric and made to volume with ACN

≠

500 mL H₂O added to 500 mL ACN



HPLC System: Agilent 1100 with quaternary pump
Column: ZORBAX Eclipse XDB-C8 Rapid-Resolution (3.5 μ m), 4.6 x 50 mm
Agilent Part No. 935967-906
Mobile Phases: Dial-a-Mix= A: water B: MeOH, pump 50% B
Vol/Vol=250 mL water + 250 mL MeOH, pump 100%
Mix-to-Mark = 250 mL MeOH, fill to 500 mL with water, pump 100%
Premixed (w/w) = 200 g MeOH + 200 g water, pump 100%
Detection: UV 254 nm
Flow: 1 mL/ min.
Temperature: ambient



1. Uracil
2. Butylparaben
3. Napthalene
4. Dipropylphthalate
5. Acenaphthene

- Method used to prepare mobile phase can significantly affect the elution
- **Be consistent and document the process**



Effect of Mobile Phase Preparation on Chromatography, 5988-6476EN

What Should I Test to Make a Robust Method?

pH 4.5 shows selectivity change from lot-to-lot for basic compounds



pH 3.0 shows no selectivity change from lot-to-lot



For method ruggedness


- Test three different column lots
- Compare R_s for the three lots
 - If ΔR_s is too large, modify method

Store RPLC Columns in 100% Acetonitrile When Not in Use

Storage solvent:


95/5 H₂O/CH₃CN

+ 0.1% **Formic Acid**




45/55 H₂O/CH₃CN

(no modifier)



100% CH₃CN

76% 0.1% FA in H₂O, 24% CH₃CN, 0.4 mL/min, isocratic, 2.1 x 150 mm columns, 60 °C, DAD: 254 nm, 80 Hz, Sample: uracil, maleic acid, imipramine, amitriptyline, methyl paraben, acetophenone




Method Setup

- What method parameters should I optimize?
- Should I use default values?



Instrument Configuration



Dwell Volume: from formation of gradient to top of column

-minimize for faster equilibration and more efficient gradient formation


Extracolumn volume from injection to detector (flow cell) outside of the column

Minimize to reduce band broadening, for sharper peaks and better resolution

Comparison of Gradient Delay Volume (Dwell Volume)


1290 Infinity II Flexible Pump (Quaternary)

- Integrated degasser
- Four solvent channels with concurrent mixing of all four channels
- Lower in price, typically, than binary pump




1290 Infinity II High Speed Pump (Binary)


- Integrated degasser
- Four solvent channels available, mixing of two channels possible
- Better performance concept is widely accepted
- Greater control over dwell volume compared to quaternary pump



Delay Volume Profiles



Chromatographic Test Results with Different Delay Volumes



Dispersion Reduces HPLC Performance

What is dispersion?

- Original sample concentration being diluted as it is carried through the system plumbing (extracolumn volume)


What increases dispersion?

- Connecting tubing that is too long
- Connecting tubing that is too large in diameter
- Connections that have gaps and form small mixing chambers



$$\sigma_{v,\text{ext}}^2 = \frac{\pi d^4 L_{cap} u_{cap}}{96D_m}$$

$\sigma_{v,\text{ext}}^2$ is the volume variance
d is the tubing diameter
L is the tubing length
u is the linear velocity of the liquid
 D_m is the molecular diffusion coefficient

Optimizing Connecting Tubing Volume for UHPLC Columns



Length	10 mm	50 mm	100 mm	150 mm
Tubing id	Volume	Volume	Volume	Volume
0.17 mm (green)	0.227 μ L	1.1 μ L	2.27 μ L	3.3 μ L
0.12 mm (red)	0.113 μ L	0.55 μ L	1.13 μ L	1.65 μ L



Pump Setting

Method of G7104A (DEBA300770) Quat. Pump (G7104A)

Flow
1.000 mL/min

Solvents
 Enable Blend Assist

A: 90.00 % 100.0 % Water V.03
B: 10.00 % 100.0 % Acetonitrile V.03
C: 0.00 % 100.0 % Acetonitrile V.03
D: 0.00 % 100.0 % Water V.03

Pressure Limits
Min: 0.00 bar Max: 1,300.00 bar

Stoptime **Posttime**
 As Injector/No Limit Off
 3.00 min 1.50 min

Advanced

Minimum Stroke
 Automatic
 20.00 µL

Compressibility
 Use Solvent Types Slow down for pressure sensitive columns

Maximum Flow Gradient
Flow ramp up: 100.000 mL/min² Flow ramp down: 100.000 mL/min²

Primary Channel
Automatic

Mixer Selection
Use Mixer if installed

▶ Timetable (1/100 events)
▶ ISET

Ok Apply Cancel

Optimize Autosampler Performance

Reduce sample carryover

Method of G4226A (DE93000256)

Injection

Injection volume: 2.00 μ l

Needle wash

Enable Needle Wash

Mode: Flush Port

Time: 50 s

Location:

Repeat: 3

Stoptime Posttime

Improved accuracy for chilled samples

Draw speed: 100.0 μ l/min (Default 200ul/min)

Eject speed: 200.0 μ l/min

Draw position: 0.0 mm

Equilibration time: 1.2 sec

Sample flush out factor: 5.0 times injection volume

Vial/Well bottom sensing

Automatic delay volume reduction

Enable overlapped injection

When Sample Is Flushed Out

After Period Of Time

0.00 min

Optimize Autosampler Performance – Draw Position/Bottom Sensing

Needle Height Position

Offset: mm

Use Vial/Well Bottom Sensing



Draw position: mm

Equilibration time: sec

Sample flush out factor: times injection volume

Vial/Well bottom sensing

Draw Position/Needle Height Position Offset = 0	Vialsampler G1329B/G7129A/B	Well-plate sampler G1367E/G4226A	Multisampler G7167A/B
	2 mL vial (sample tray)	2 mL vial 54 vial tray	2 mL vial 54 vial tray
Without bottom sensing	2 mm	4 mm	5 mm
With bottom sensing	x	1 mm	2 mm



VWD and DAD Settings

Wavelength: 250 nm
Peakwidth: > 0.1 min (2 s resp. time) (5 Hz)

Zero Offset: 5 %
Attenuation: 1000 mAU

Signal Polarity
Positive (+)
Negative (-)

Miscellaneous
Lamp on required for acquisition
Scan Range: 190 to 200 nm
Step: 2 nm

Additional Signals
Acquire Signal without Reference
Acquire Reference only

No bandwidth setting
No slit width setting

Only use reference or not option

Signal	Acquire	Wavelength	Bandwidth	Reference Wavelength	Reference Bandwidth
Signal A	<input checked="" type="checkbox"/>	254.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm
Signal B	<input checked="" type="checkbox"/>	254.0	4.0	<input type="checkbox"/> 360.0	100.0 nm
Signal C	<input type="checkbox"/>	214.0	4.0	<input type="checkbox"/> 360.0	100.0 nm
Signal D	<input type="checkbox"/>	230.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm
Signal E	<input type="checkbox"/>	260.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm
Signal F	<input type="checkbox"/>	273.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm
Signal G	<input type="checkbox"/>	280.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm
Signal H	<input type="checkbox"/>	250.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0 nm

Peakwidth
> 0.013 min (0.25 s response time) (20 Hz)
< 0.0008 min (0.008 s response time) (240 Hz)
> 0.0008 min (0.016 s response time) (240 Hz)
> 0.0016 min (0.031 s response time) (160 Hz)
> 0.0031 min (0.063 s response time) (80 Hz)
> 0.0063 min (0.13 s response time) (40 Hz)
> 0.013 min (0.25 s response time) (20 Hz)
> 0.025 min (0.5 s response time) (10 Hz)
> 0.05 min (1 s response time) (5 Hz)
> 0.1 min (2 s response time) (2.5 Hz)
> 0.2 min (4 s response time) (1.25 Hz)
> 0.4 min (8 s response time) (0.62 Hz)
> 0.85 min (16 s response time) (0.31 Hz)

Store: All
Range from: 190.0 to 400.0 nm
Step: 2.0 nm

Zero Offset: 5 %
Attenuation: 1000 mAU

Margin for negative Absorbance: 100 mAU
Slit: 4 nm

Autobalance
Prerun
Postrun



Lamps on required for acquisition
UV Lamp

T timetable (empty)

Ok Apply Cancel



DAD Setting – Choose the Right Signal and Reference

Signals					
	Acquire	Wavelength	Bandwidth	Reference Wavelength	Reference Bandwidth
Signal A	<input checked="" type="checkbox"/>	254.0	4.0	<input checked="" type="checkbox"/>	360.0
Signal B	<input checked="" type="checkbox"/>	254.0	4.0	<input type="checkbox"/>	100.0



Gradient: 10–100% ACN in 3 min




DAD Setting – Choose the Right Signal and Reference



Reference wavelength may not be necessary


DAD Setting – Choose the Right Signal and Reference

Move away from the UV cutoffs of mobile phase/additive




What Data Rate Should I Choose?

InfinityLab Poroshell 120 EC-C18, 2.1 × 50 mm, 1.9 μm
20 mM sodium phosphate pH 7 in water with acetonitrile premixed 40/60
0.5 mL/min




App note: 5991-7560EN



Fast data collection rates must be used with Agilent InfinityLab Poroshell 1.9 μm columns to accurately measure the efficiency of the column, especially for early eluting compounds such as butylparaben ($k' = 1.3$).

DAD Setting – Choose the Right Sampling Rate



Do not use peak width smaller than necessary

Peakwidth

Stoptime

As P min

> 0.013 min (0.25 s response time) (20 Hz)
< 0.0008 min (0.008 s response time) (240 Hz)
> 0.0008 min (0.016 s response time) (240 Hz)
> 0.0016 min (0.031 s response time) (160 Hz)
> 0.0031 min (0.063 s response time) (80 Hz)
> 0.0063 min (0.13 s response time) (40 Hz)
> 0.013 min (0.25 s response time) (20 Hz)
> 0.025 min (0.5 s response time) (10 Hz)
> 0.05 min (1 s response time) (5 Hz)
> 0.1 min (2 s response time) (2.5 Hz)
> 0.2 min (4 s response time) (1.25 Hz)
> 0.4 min (8 s response time) (0.62 Hz)
> 0.85 min (16 s response time) (0.31 Hz)

Column: ZORBAX Eclipse Plus C18, 2.1 x 50 mm, 1.8 μ m
Column temperature: 35 $^{\circ}$ C; flow rate: 1 mL/min
Gradient: 10-100% ACN in 3 min
Signal: 254 nm, bandwidth: 4 nm
Reference: 360 nm, bandwidth: 100 nm

Agilent InfinityLab

Designed to seamlessly integrate into the InfinityLab family

You can rely on **Agilent InfinityLab** LC instruments, columns, and supplies to deliver rugged quality and robust analytical results.

www.agilent.com/chem/infinitylab

InfinityLab columns



InfinityLab supplies



InfinityLab
LC series


InfinityLab
LC/MSD series




InfinityLab
Accessories

Agilent Resources for Support

- Resource page <http://www.agilent.com/chem/agilentresources>
 - Quick reference guides, product catalogs
 - Online selection tools, “How-to” videos
 - Column user guides - <https://www.agilent.com/en-us/support/liquid-chromatography/kb005965>
 - Biocolumn user guides - <https://www.agilent.com/en/support/liquid-chromatography/kb005960>
- Tech support: <http://www.agilent.com/chem/techsupport>
- InfinityLab LC Supplies catalog ([5991-8031EN](#))
- Agilent University <http://www.agilent.com/crosslab/university>
- YouTube – [Agilent Channel](#)
- Your local product specialists
- Subscribe to Agilent Peak Tales podcasts at peaktales.libsyn.com



Contact Agilent Chemistries and Supplies Technical Support



Available in the USA and Canada 8-5 all time zones

1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com