

Mass spectrometry based assessment chimeric mouse liver metabolite profiles following oral dosing of troglitazone

IMSC 2012 PWe-093

Alan Barnes¹; Neil J Loftus¹; Kirsten Hobby¹;
Ian Wilson²; Yoshio Morikawa³

¹Shimadzu MS/BU, Manchester, UK;

²Astra Zeneca, Alderley Park, Cheshire, UK;

³PhoenixBio Co. Ltd, Higashi-Hiroshima, Japan

Mass spectrometry based assessment chimeric mouse liver metabolite profiles following oral dosing of troglitazone

Introduction

In the continuing search for new chemical entities the use of chimeric mice with humanized livers are being used in the search for unexpected drug metabolites. Chimeric mice, in which the majority of the hepatocyte population of the mouse liver has been replaced by human hepatocytes, have the capacity to express human Phase I and II metabolic enzymes and hepatic transporter proteins with gene-expression profiles and phenotypes similar (up to 85%) to those of the original donor liver. To assess the viability of the chimeric Phoenix Bio (PXB) mouse in

modeling human liver metabolism, troglitazone (TGZ) was dosed orally over 7 days at two dose concentrations (300 & 600 mg/kg). In pre-clinical studies TGZ showed inter-species differences in metabolism particularly in sulfation and glucuronidation pathways. The present study evaluated the metabolic profile of troglitazone and endogenous metabolites in the PXB compared to control mice (severe combined immunodeficiency - SCID) using high mass accuracy MS/MS analysis.

Materials and Methods

Liver extracts from SCID (control) and PXB (chimeric) mice were analyzed using a high resolution LC/MSn system (Nexera LC coupled with a LCMS-IT-TOF; Shimadzu Corporation). Both aqueous and organic extracts were analyzed using a Phenomenex Kinetex column (C18 1.7 μ m, 2.1 \times 100 mm); aqueous components were separated at a flow rate of 0.6 mL/min, with the column maintained at 30°C. The chromatographic system used a binary solvent system delivered as a gradient of solvent A (H₂O containing 0.1% formic acid, 10mM ammonium acetate) and solvent B (ACN containing 0.1% formic acid). The gradient conditions were: 5% B (5 min), to 35% (3 min), to 50% (22 min), to 95% (2.5 min) held for 5 min, re-equilibration 2.5 min. The solvent composition was then held at 100% B for 6.5 min after which the column was returned to 5% B

over the next 2.5 min, making a total cycle time of 40 min per sample; organic extracts were separated using a different gradient method (Castro-Perez *et al.* 2010) – (mobile phase: A – water:acetonitrile (60:40) 10 mM ammonium formate pH5, B - propan-2-ol:acetonitrile (90:10) 10 mM ammonium formate) at a flow rate of 0.5 mL/min. The LCMS-IT-TOF acquired positive and negative MS and MS² data using high speed polarity switching (*m/z* 150-1250). Profiling Solution software (Shimadzu, Japan) was applied to metabolite profiling analysis to assess the impact of changes in lipid profiles. MetID Solution software was used for a targeted metabolomics study in addition to searching for drug metabolites through similarity scoring MSⁿ data of potential metabolites by comparing common fragment ions to parent drug MSⁿ data.




Fig. 1 Liver metabolite profile of control mouse (SCID) following oral administration of troglitazone (TGZ).


Mass spectrometry based assessment chimeric mouse liver metabolite profiles following oral dosing of troglitazone

Results

Troglitazone metabolism




Analysis of aqueous liver extracts by accurate mass negative ion MSⁿ enabled detection of metabolites by MetID Solution software (Fig. 1). Confirmation of troglitazone

metabolites was also possible through analysis of common fragmentation data (Fig. 2).





| MS/MS fragments | Fragment Formula | Expected mass | Mass accuracy (ppm) |
|-----------------|---|---------------|---------------------|
| 412.1581 | C ₂₃ H ₂₇ NO ₄ S | 412.1588 | -1.7 |
| 397.1488 | C ₂₃ H ₂₆ O ₄ S | 397.1479 | 2.3 |
| 369.1546 | C ₂₂ H ₂₆ O ₃ S | 369.1530 | 4.3 |
| 363.1577 | C ₂₃ H ₂₄ O ₄ | 363.1602 | -6.9 |
| 276.0714 | C ₁₄ H ₁₅ NO ₃ S | 276.0700 | 5.1 |
| 179.0175 | C ₉ H ₈ O ₂ S | 179.0172 | 1.7 |
| 162.0576 | C ₉ H ₉ NO ₂ | 162.0561 | 9.3 |
| 145.0299 | C ₉ H ₆ O ₂ | 145.0295 | 2.8 |

MS/MS Troglitazone: C₂₄H₂₇NO₅S [M-H] 440.1537 RT 21.79 min



MS/MS TGZ-sulfate: C₂₄H₂₇NO₈S₂ [M-H] 520.1105 RT 13.37 min



MS/MS TGZ- glucuronide: C₃₀H₃₅NO₁₁S [M-H] 616.1858 RT 11.07 min

Fig. 2 Fragmentation analysis of troglitazone by accurate mass MSⁿ data.

Common fragment ions and neutral loss information consistent to troglitazone parent enabled characterization of metabolite structures.

Mass spectrometry based assessment chimeric mouse liver metabolite profiles following oral dosing of troglitazone

Table 1 Averaged peak area data of troglitazone and metabolites detected in aqueous liver extracts.

| Peak ID | Assignment | MS2 | RT | <i>m/z</i> [M-H] ⁻ | SCID 600 mg | PXB 600 mg |
|--------------|------------------------|-----|-------|-------------------------------|-------------|------------|
| Troglitazone | Parent | + | 21.79 | 440.1537 | 2,807,994 | 3,898,820 |
| M1 | Di-hydroxy glucuronide | | 8.61 | 648.1756 | 72,254 | 71,725 |
| M2 | Hydrated glucuronide | + | 8.40 | 634.1963 | 4,842,608 | 3,460,630 |
| M3 | Hydrated sulfate | + | 8.81 | 538.1211 | 5,390,498 | 6,246,988 |
| M4 | Hydroxy sulfate | + | 9.18 | 536.1054 | 25,491 | 28,861 |
| M9 | Di-hydroxy | + | 9.63 | 472.1435 | 177,029 | 184,855 |
| M10 | Hydroxy glucuronide | + | 9.88 | 632.1807 | 174,260 | 125,705 |
| M12 | Hydroxy sulfate | + | 11.59 | 536.1054 | 336,482 | 230,705 |
| M13 | Glucuronide | + | 11.07 | 616.1858 | 12,618,486 | 8,414,646 |
| M15 | Sulfate | + | 13.37 | 520.1105 | 25,852,882 | 26,671,871 |
| M16 | Di-hydroxy | + | 10.52 | 472.1435 | 345,952 | 150,980 |
| M18 | Di-hydroxy | | 11.12 | 472.1435 | 105,877 | 75,805 |
| M27 | Mono-hydroxy | + | 14.77 | 456.1486 | 193,239 | 280,266 |
| M30 | Mono-hydroxy | + | 15.74 | 456.1486 | 2,657,528 | 2,307,728 |

Peak area data comparing relative levels of troglitazone metabolites showed differences in metabolic profiles were also observed between PXB and SCID mice; consistent with metabolic profiles reported in human and mouse, the sulfate conjugate was more abundant in PXB than SCID whilst the glucuronide metabolite was greater in mouse.

Endogenous metabolite profiling

Organic liver extracts were analyzed to examine endogenous lipid differences between PXB and SCID livers. Data was aligned using Profiling Solution software (Shimadzu Corporation) and principal component analysis (PCA) was performed to examine group differences using Simca-P (Umetrics).




Fig. 3 Statistical analysis of organic liver extracts comparing all SCID to all PXB samples.
a) PCA analysis revealed two main experimental groups (PXB and SCID) with no clear grouping associated with dosing of troglitazone. Tight clustering of QA/QC samples indicated good system stability throughout the sample analysis period.

b) OPLS-DA S-plot analysis comparing PXB to SCID enabled ions of highest significance to be identified. Two diacylglycerophosphocholine compounds (labeled) were detected at significantly higher levels in PXB mice compared to SCID

Mass spectrometry based assessment chimeric mouse liver metabolite profiles following troglitazone

MetID Solution was used to perform a targeted search of known endogenous metabolites using LipidMaps entry information from the following compound classes: phosphatidic acid, phosphatidylglycerol, phosphatidylserine, phosphatidylethanolamine, phosphatidylinositol, phosphatidylcholine. The analysis enabled identification of over 80 ions that differed significantly between sample groups (concise summary: Table 2). Putative identifications were made based on mass accuracy and isotope score. Fold differences are shown between PXB and SCID at no dose (0 mg/kg), high dose (600 mg/kg) and for all animals averaged (0, 300 and 600 mg/kg). Although the aim of the data

analysis was to identify compounds that differed between PXB and SCID mice, the data analysis also revealed subtle differences occurring possibly as a result of troglitazone dosing. Some compounds such as the glycerophosphocholine compounds were consistent in up or down regulation irrespective of dosing, hence showing most significance in S-plot analysis (Fig. 3b) due to homogenous variance averaged across all dosing groups. Conversely other lipid species exhibited differences in the fold change, although still up or down due to being PXB or SCID, show that administration of troglitazone may influence the concentration.

Table 2 Endogenous metabolites identified as significantly increased (green) or decreased (red) in PXB mice compared to SCID mice at 0 mg, 600 mg dosing and data from all animals averaged (SCID- indicates not detected in SCID)

| DB reference | Putative ID | Formula | Ion | m/z | RT | % RSD | 0 mg | | 600 mg | | All animals | |
|---------------|---|-------------|---------------------|----------|-------|-------|--------|------------|------------|------------|-------------|------------|
| | | | | | | | QA/QC | PXB / SCID | PXB / SCID | PXB / SCID | PXB / SCID | PXB / SCID |
| LM GL03010065 | TG(16:0/16:1(9Z)/18:3(9Z,12Z,15Z)) | C53H94O6 | [M+H] ⁺ | 827.7123 | 24.35 | 8.2 | 5.72 | SCID- | SCID- | | | |
| LM GP01010395 | PC(10:0/20:0) | C38H76NO8P | [M+H] ⁺ | 706.5381 | 9.93 | 6.1 | 14.61 | 10.62 | 10.94 | | | |
| M ID370 | Glycerophosphocholine | C8H20NO6P | [M+H] ⁺ | 258.1101 | 0.48 | 3.8 | 22.48 | 8.08 | 10.78 | | | |
| LM GL03010078 | TG(16:1(9Z)/16:1(9Z)/18:3(9Z,12Z,15Z)) | C53H92O6 | [M+H] ⁺ | 825.6967 | 24.16 | 7.0 | 7.92 | 20.67 | 10.46 | | | |
| LM GL03010018 | TG(16:1(9Z)/14:0/18:1(9Z)) | C51H84O6 | [M+NH] ⁺ | 820.7389 | 24.12 | 4.8 | 11.02 | 8.59 | 8.76 | | | |
| LM GP01010490 | PC(14:0/18:1(11Z)) | C40H78NO8P | [M+H] ⁺ | 732.5558 | 11.21 | 13.3 | 20.25 | 9.47 | 8.76 | | | |
| HM DB01235 | 5-Aminoimidazole ribonucleotide | C8H14N3O7P | [M+H] ⁺ | 296.0642 | 0.41 | 3.2 | 11.48 | 7.53 | 8.26 | | | |
| LM GP10020005 | PA(0-16:0/14:1(9Z)) | C33H65O7P | [M+H] ⁺ | 605.4541 | 14.08 | 6.9 | 10.93 | 6.63 | 6.65 | | | |
| LM GP10010088 | PA(13:0/22:2(13Z,16Z)) | C38H71O8P | [M-H] ⁻ | 685.4814 | 10.26 | 7.2 | 15.10 | 2.98 | 6.10 | | | |
| LM GP06010075 | PI(14:0/22:2(13Z,16Z)) | C45H83O13P | [M-H] ⁻ | 861.5499 | 10.89 | 4.7 | 4.81 | 4.64 | 5.82 | | | |
| LM ST05040015 | Tauroursodeoxycholic acid | C26H45NO6S | [M-H] ⁻ | 498.2895 | 6.60 | 3.2 | 9.79 | 4.34 | 5.41 | | | |
| LM GP04020069 | PG(0-20:0/22:0) | C48H97O9P | [M+H] ⁺ | 849.6943 | 24.00 | 7.1 | 4.83 | 5.57 | 5.37 | | | |
| LM GP10020004 | PA(0-16:0/14:0) | C33H67O7P | [M+H] ⁺ | 607.4697 | 17.11 | 6.8 | 6.24 | 5.72 | 5.36 | | | |
| LM GP01010508 | PC(14:0/20:5(5Z,8Z,11Z,14Z,17Z)) | C42H74NO8P | [M+H] ⁺ | 752.5225 | 8.50 | 8.2 | 6.73 | 4.40 | 4.59 | | | |
| LM GL03010166 | TG(17:2(9Z,12Z)/17:2(9Z,12Z)/18:2(9Z,12Z)) | C55H94O6 | [M+H] ⁺ | 851.7123 | 24.18 | 5.7 | 5.91 | 4.36 | 4.56 | | | |
| LM GP01010490 | PC(14:0/18:1(11Z)) | C40H78NO8P | [M+H] ⁺ | 732.5538 | 10.23 | 4.1 | 6.65 | 4.53 | 4.51 | | | |
| LM GP01010512 | PC_LM GP01010512 | C44H76NO8P | [M+H] ⁺ | 778.5381 | 8.62 | 9.0 | 6.48 | 5.67 | 4.45 | | | |
| LM GP01010490 | PC(14:0/18:1(11Z)) | C40H78NO8P | [M+H] ⁺ | 732.5538 | 10.19 | 3.8 | 6.41 | 4.41 | 4.38 | | | |
| LM GL03010140 | LM GL03010140 | C55H96O6 | [M+H] ⁺ | 853.7280 | 24.36 | 4.5 | 4.59 | 4.68 | 4.23 | | | |
| LM GP01010494 | PC(14:0/18:2(11Z,14Z)) | C40H76NO8P | [M+H] ⁺ | 730.5381 | 9.21 | 3.7 | 5.55 | 4.39 | 4.18 | | | |
| LM GP01010541 | PC(15:0/18:1(11Z)) | C41H80NO8P | [M+H] ⁺ | 746.5694 | 15.46 | 4.1 | 4.23 | 4.27 | 4.06 | | | |
| LM GP01010633 | PC(16:0/20:5(5Z,8Z,11Z,14Z,17Z)) | C44H78NO8P | [M+H] ⁺ | 780.5538 | 9.51 | 3.2 | 2.57 | 2.54 | 2.33 | | | |
| LM GP01050125 | PC(16:1(9Z)/0:0) | C23H46NO7P | [M+H] ⁺ | 480.3085 | 5.68 | 4.4 | -3.63 | -1.80 | -2.30 | | | |
| LM GP01010645 | PC(16:0/22:5(4Z,7Z,10Z,13Z,16Z)) | C46H82NO8P | [M+H] ⁺ | 808.5851 | 10.50 | 14.0 | -3.39 | -1.88 | -2.34 | | | |
| LM GP06010076 | PI(14:0/22:4(7Z,10Z,13Z,16Z)) | C45H79O13P | [M-H] ⁻ | 857.5166 | 9.36 | 5.0 | -2.62 | -2.24 | -2.38 | | | |
| LM GP10020032 | PA(0-18:0/18:4(6Z,9Z,12Z,15Z)) | C39H71O7P | [M+H] ⁺ | 683.5010 | 18.30 | 7.6 | -3.24 | -2.20 | -2.46 | | | |
| LM GL02010197 | DG(20:3(8Z,11Z,14Z)/20:4(5Z,8Z,11Z,14Z)/0:0) | C43H70O5 | [M+H] ⁺ | 667.5296 | 18.13 | 2.9 | -2.95 | -2.43 | -2.49 | | | |
| LM GL03010722 | TG(18:3(9Z,12Z,15Z)/18:3(9Z,12Z,15Z)/20:1(11Z)) | C59H100O6 | [M+H] ⁺ | 905.7593 | 18.37 | 7.0 | -10.88 | -1.10 | -3.98 | | | |
| LM GP01011755 | PC(19:0/22:6(4Z,7Z,10Z,13Z,16Z,19Z)) | C49H86NO8P | [M+H] ⁺ | 848.6164 | 13.07 | 3.5 | -6.06 | -3.26 | -4.47 | | | |
| LM GP01011670 | PC(18:3(6Z,9Z,12Z)/22:6(4Z,7Z,10Z,13Z,16Z,19Z)) | C48H78NO8P | [M+H] ⁺ | 828.5538 | 9.53 | 9.3 | -5.62 | -3.80 | -4.64 | | | |
| LM GP01010788 | PC(18:0/20:2(11Z,14Z)) | C46H88NO8P | [M+H] ⁺ | 814.6320 | 15.13 | 4.0 | -5.38 | -3.58 | -4.71 | | | |
| LM GP01020004 | PC(0-10/18:0) | C25H52NO7P | [M-H] ⁻ | 508.3409 | 6.35 | 17.7 | -10.58 | -3.88 | -4.85 | | | |
| LM GP01010788 | PC(18:0/20:2(11Z,14Z)) | C46H88NO8P | [M+H] ⁺ | 814.6320 | 15.54 | 5.7 | -26.45 | -4.61 | -5.38 | | | |
| LM GP03010199 | PS(16:0/22:1(11Z)) | C44H84NO10P | [M+H] ⁺ | 818.5906 | 8.51 | 6.2 | -50.23 | -5.10 | -5.65 | | | |
| LM GP01011028 | PC(20:0/22:6(4Z,7Z,10Z,13Z,16Z,19Z)) | C50H98NO8P | [M+H] ⁺ | 862.6320 | 14.48 | 4.7 | -20.00 | -9.37 | -14.21 | | | |
| LM GP01020072 | PC(0-16:0/4:0) | C28H58NO7P | [M+H] ⁺ | 552.4024 | 7.10 | 6.2 | -43.04 | -12.65 | -14.76 | | | |

Conclusion

- Human specific troglitazone metabolism, consistent to published data, was shown from PXB mice.
- Endogenous lipid differences between PXB and SCID were detected, some consistent irrespective of troglitazone dosing and others that may be influenced by troglitazone dosing.
- MetID Solution combined use of accurate mass and isotope scoring enabled greater confidence in putative metabolite identification.



For Research Use Only. Not for use in diagnostic procedures.
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.