



# Characterizing LNP mRNA

LC tools for purity, identity, integrity, and concentration determining measurements




# A Changing Industry

mRNA is becoming a highly effective drug substance for both *in vivo* and *ex vivo* treatment of human infections, diseases and cancers. When encapsulated into a lipid nanoparticle, mRNA can be efficiently delivered to patient cells to be expressed into vaccine antigens, enzyme replacements to bolster a patient's own gene expression, or someday a gene editing apparatus to alter a patient's somatic genome. The FDA approved COVID-19 mRNA vaccines have provided a successful formula for creating potent LNPs. With today's techniques, four types of lipids encase the mRNA. A synthetic ionizable lipid is combined with a zwitterionic phosphatidylcholine lipid, cholesterol and a PEGylated lipid to yield 700–800 Å diameter nanoparticles. mRNA payloads are incorporated that range from 1,000 to up 10,000 nucleotides in length. Given their complexity, these LNP RNA drug products must be comprehensively characterized and tested. Chromatography can expedite the needed assays and thereby confirm physicochemical properties of the intact particle, check their compositions and raw materials for impurities, and confirm the molecular integrity and modifications of the mRNA drug substance.



Structure of lipid Nanoparticle



## Chromatography and Method Options

Size Exclusion

Anion Exchange

Poly A Tail LC

Oligonucleotide Mapping

5' Cap Analysis

Lipid Testing

## Waters Peer Reviewed Articles

Waters scientists and collaborators are publishing on this subject. Make sure to visit the Resource Tab on our [waters.com/GTx](https://www.waters.com/GTx) website to keep up to date on the literature.

## SEC to Measure Integrity of mRNA Drug Substance

Size Exclusion Chromatography (SEC) is a widely employed separation technique for isolating species based upon differences in hydrodynamic radius. SEC of larger nucleic acids requires novel method development to test and report on critical quality attributes (CQA) as they relate to the safety and efficacy of the drug. A GTxResolve™ Premier BEH™ SEC 450 Å 2.5 µm Column provides high resolution separations of small to medium sized nucleic acids species in various mobile phase conditions. New SEC separation capabilities with denaturing mobile phases have also come to be established. Formulated LNP samples can be injected onto an SEC column running a buffer comprised of 0.2% SDS + 20% isopropanol to facilitate quick measurements on the nucleic acid payloads within (Figure 2).




Figure 1. Components of 6000 nt RNA ladder were resolved using a Waters GTxResolve Premier BEH 450 Å SEC Column as reported in Waters Application Note: [720008061](#).




Figure 2. Denaturing SEC performed with a GTxResolve Premier BEH SEC 450 Å 2.5 µm Column and a mobile phase containing 0.2% SDS and 20% isopropanol. Injection of Comirnaty™ COVID-19 mRNA Vaccine 10 µg/dose drug product with simultaneous on-column dissolution of the LNP and detection of the mRNA payload.

## Measuring Intact LNP and Large Nucleic Heterogeneity by SEC

Size exclusion chromatography (SEC) is also becoming a powerful technique to assess the heterogeneity and integrity of large nucleic acids and their corresponding LNPs. Contaminating impurities such as truncated and aggregated mRNA have implications on safety and efficacy and must therefore be monitored. In addition, first-of-their kind measurements on the heterogeneity of LNPs are now becoming possible such that new structure-function relationships can be elucidated. GTxResolve Premier SEC 1000 Å 3 µm Columns are built with MaxPeak Premier High Performance Surface Hardware and a novel particle technology to ensure that efficient, high recovery separations can be readily achieved and that new CQAs can be quickly reported.




Figure 3. GTxResolve Premier SEC 1000 Å 3 µm Columns contain high efficiency particles modified with a novel ethylene bridged HO-PEO surface chemistry, which produces high recovery analyses of 200 to 1000 Å diameter GTx drugs, such as mRNA, LNPs, and viral vectors. This crosslink protected, hydrophilic surface chemistry also ensures long column lifetimes, low secondary interactions and improved MALS sensitivity.




Figure 4. Comirnaty COVID-19 mRNA Vaccine 10 µg/dose drug product and intact LNP heterogeneity measurement using a GTxResolve Premier SEC 1000 Å 3 µm 4.6 x 150 mm Column.





## 5' Cap Analysis and Oligo Mapping Analysis

An mRNA's sequence and its modifications need to be confirmed to ensure the mRNA will reach its target efficacy. The 5' cap can be characterized after sample prep with DNA probes and RNase H which is specific to DNA/RNA duplexes. Further sequence confirmation can be obtained with mRNA digestion leveraging residue specific endonucleases. Oligo batch tested BEH C<sub>18</sub> sorbents provide rugged columns for each of the downstream separations. MaxPeak™ High Performance Surfaces ensure quick method starts by eliminating conditioning effects and improving analyte recovery.




Figure 6. 5' cap digestion product from RNase H cleavage (A) and RNase T1 oligonucleotide map of a tail-less luciferase mRNA (B) as published in Waters Application Notes: [720007329](#) and [720007669](#).

# Poly A Tail Analysis by RPLC

The length and structure of a 3' poly A tail must be optimized to confer desired half life and ribosome binding affinity properties to an mRNA. RNase cleavage can be applied to digest an mRNA down to its poly A motif. An oligonucleotide batch tested widepore BEH 300 Å C<sub>18</sub> column provides a high resolution separation of the liberated poly A tail, and the use of strong ion pairing agents makes it possible to achieve single residue resolution.




Figure 7. IP-RPLC chromatogram of a T1 digested mRNA sample and its constituent poly A tail as obtained with an ACQUITY Premier Oligonucleotide BEH 300Å 1.7 µm column (Waters Application Note: [720007873](#); [Anal. Chem. 2023, 95, 38, 14308–14316](#)).

# Ultrafast RPLC to Increase Throughput

Your assay could achieve a new level of throughput and productivity with adoption of short bed length columns. Waters has developed low adsorption ACQUITY™ Premier 2.1 x 20 mm Columns packed with 1.7 µm oligonucleotide batch tested and selected stationary phases. These columns have been carefully optimized for both efficiency and mechanical durability, offering improved throughput for purity, identity and critical quality attribute testing. A key insight from recently published work shows that throughput can be significantly improved by selecting the appropriate column size. The larger the oligonucleotide or mRNA digestion component the more likely they are to exhibit a bind and elute chromatography mechanism.




Figure 8. Isotherms of a 20 and 100-mer oligonucleotide that demonstrate the quick change of the analyte from an adsorbed to desorbed state upon an increase in eluent strength. A example 1-minute IP-RPLC analysis. ([Anal. Chem. 2023, 95, 27, 10448–10456](#)).



# Ordering Information

The columns, standards, and reagents that can help you characterize LNP mRNA are provided below.



*Quality and Design Matters:  
Use consumables built and quality control tested for the highly ionic analytes comprised within cell and gene therapy drug substances and products.*



*MaxPeak High Performance Surfaces are based on a vapor deposited organosilica material that is incorporated into Waters column hardware to yield higher analyte recoveries and improved peak shapes.*

## SEC

### GTxResolve Premier BEH SEC 450 Å 2.5 µm Columns

|                 | Dimension | 150 mm                    | 300 mm                    | Guard                     |
|-----------------|-----------|---------------------------|---------------------------|---------------------------|
| Standard Column | 4.6 mm    | <a href="#">186010584</a> | <a href="#">186010585</a> | <a href="#">186010583</a> |
|                 | 7.8 mm    | <a href="#">186010586</a> | <a href="#">186010587</a> | <a href="#">186010583</a> |

### GTxResolve Premier SEC 1000 Å 3 µm Columns

|                 | Dimension | 150 mm    | 300 mm    | Guard     |
|-----------------|-----------|-----------|-----------|-----------|
| Standard Column | 4.6 mm    | 186010735 | 186010736 | 186010733 |
|                 | 7.8 mm    | 186010737 | 186010738 | 186010733 |



## AEX

### Anion Exchange Columns

|                            | Dimension    | P/N                       |
|----------------------------|--------------|---------------------------|
| Protein-Pak HiRes Q Column | 4.6 x 100 mm | <a href="#">186004931</a> |
| Gen-Pak FAX Column         | 4.6 x 100 mm | <a href="#">WAT015490</a> |



## STANDARDS AND REAGENTS

### Standards and Reagents

|                             | P/N                       |
|-----------------------------|---------------------------|
| dsDNA 50 to 1350 Ladder     | <a href="#">186010778</a> |
| ssDNA 10 to 60 Ladder       | <a href="#">186009449</a> |
| ssDNA 20 to 100 Ladder      | <a href="#">186009448</a> |
| ssDNA 20-mer LC-MS Standard | <a href="#">186009451</a> |





## OLIGO RPLC

ACQUITY Premier Oligonucleotide BEH C<sub>18</sub> 1.7 µm Columns

|                      | Diameter | 130 Å                     |                           |                           | 300 Å                     |                           |                           |
|----------------------|----------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
|                      |          | 50 mm                     | 100 mm                    | 150 mm                    | 50 mm                     | 100 mm                    | 150 mm                    |
| Standard Column      | 2.1 mm   | <a href="#">186009484</a> | <a href="#">186009485</a> | <a href="#">186009486</a> | <a href="#">186010539</a> | <a href="#">186010540</a> | <a href="#">186010541</a> |
| VanGuard™ FIT Column | 2.1 mm   | <a href="#">186010685</a> | <a href="#">186010686</a> | <a href="#">186010687</a> | <a href="#">186010754</a> | <a href="#">186010755</a> | <a href="#">186010756</a> |

XBridge™ Premier Oligonucleotide BEH C<sub>18</sub> 2.5 µm Columns

|                     | Diameter | 130 Å                     |                           |                           | 300 Å                     |                           |                           |
|---------------------|----------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
|                     |          | 50 mm                     | 100 mm                    | 150 mm                    | 50 mm                     | 100 mm                    | 150 mm                    |
| Standard Column     | 2.1 mm   | <a href="#">186009836</a> | <a href="#">186009837</a> | <a href="#">186009838</a> | <a href="#">186010542</a> | <a href="#">186010543</a> | <a href="#">186010544</a> |
|                     | 4.6 mm   | <a href="#">186009901</a> | <a href="#">186009902</a> | <a href="#">186009903</a> | <a href="#">186010545</a> | <a href="#">186010546</a> | <a href="#">186010547</a> |
| VanGuard FIT Column | 2.1 mm   | <a href="#">186010688</a> | <a href="#">186010689</a> | <a href="#">186010690</a> | <a href="#">186010757</a> | <a href="#">186010758</a> | <a href="#">186010759</a> |
|                     | 4.6 mm   | <a href="#">186010691</a> | <a href="#">186010692</a> | <a href="#">186010693</a> | <a href="#">186010760</a> | <a href="#">186010761</a> | <a href="#">186010762</a> |

ACQUITY Premier Oligonucleotide BEH C<sub>18</sub> 1.7 µm UltraFast Columns

|                 | Dimension   | 130 Å                     | 300 Å                     |
|-----------------|-------------|---------------------------|---------------------------|
| Standard Column | 2.1 x 20 mm | <a href="#">186011115</a> | <a href="#">186011021</a> |



## LIPID RPLC

## ACQUITY Premier CSH Phenyl Hexyl 1.7 µm Columns


|                     | Diameter | 130 Å                     |                           |                           |
|---------------------|----------|---------------------------|---------------------------|---------------------------|
|                     |          | 50 mm                     | 100 mm                    | 150 mm                    |
| Standard Column     | 2.1 mm   | <a href="#">186009474</a> | <a href="#">186009475</a> | <a href="#">186009476</a> |
| VanGuard FIT Column | 2.1 mm   | <a href="#">186009477</a> | <a href="#">186009478</a> | <a href="#">186009479</a> |

## LIPID RPLC


## XSelect™ Premier CSH Phenyl Hexyl 2.5 µm Columns

|                     | Diameter | 130 Å                     |                           |                           |
|---------------------|----------|---------------------------|---------------------------|---------------------------|
|                     |          | 50 mm                     | 100 mm                    | 150 mm                    |
| Standard Column     | 2.1 mm   | <a href="#">186009879</a> | <a href="#">186009880</a> | <a href="#">186009881</a> |
|                     | 4.6 mm   | <a href="#">186009886</a> | <a href="#">186009887</a> | <a href="#">186009888</a> |
| VanGuard FIT Column | 2.1 mm   | <a href="#">186009882</a> | <a href="#">186009883</a> | <a href="#">186009884</a> |
|                     | 4.6 mm   | <a href="#">186009889</a> | <a href="#">186009890</a> | <a href="#">186009891</a> |


## Application Notes




[Properties of the Gen-Pak™ FAX Column and Its Utility for Anion Exchange Analysis of Large Molecule Biologics](#)




[Synthetic mRNA Oligo-Mapping Using Ion-Pairing Liquid Chromatography and Mass Spectrometry](#)




[Salt Plug Injection Methods for Improved Anion Exchange Analyses of Large Nucleic Acids](#)




[Optimized ELSD Workflow for Improved Detection of Lipid Nanoparticle Components](#)




[Suitability of XBridge™ Premier Gx BEH™ SEC 450 Å 2.5 µm Column for Size-based Separations of Nucleic Acids](#)




[Lipid Nanoparticle Analysis: Leveraging MS to Reduce Risk](#)




[Ion-Pair Reversed-Phase Liquid Chromatography Method for Analysis of mRNA Poly\(A\) Tail Heterogeneity](#)




[Size-Exclusion Chromatography Method for Poly\(A\) Tail Analysis of mRNA](#)




[Methods for the Anion Exchange Chromatographic Analysis of mRNAs](#)



[Rapid Analysis of Lipid Nanoparticle Components Using BioAccord LC-MS System](#)




[Rapid Analysis of Synthetic mRNA Cap Structure Using Ion-Pairing RPLC with the BioAccord LC-MS System](#)




[Lipid Nanoparticle Compositional Analysis Using Charged Surface Hybrid Phenyl-Hexyl Separation With Evaporative Light Scattering Detection](#)

## Waters Journal Articles




[“Monitoring stability indicating impurities and aldehyde content in lipid nanoparticle raw material and formulated drugs”](#)

J Chrom B. 2024, 1234, 124005.




[“Challenges and emerging trends in liquid chromatography-based analyses of mRNA pharmaceuticals”](#)

J Pharm Biomed Anal. 2023, 224, 115174.




[“High-Throughput Chromatographic Separation of Oligonucleotides: A Proof of Concept Using Ultra-Short Columns”](#)

Anal. Chem. 2023, 95, 27, 10448–10456.



[“Salt gradient and ion-pair mediated anion exchange of intact messenger ribonucleic acids”](#) J Chrom Open, 2022, 2, 100031.




[“Separation of Plasmid DNA Topology Forms, Messenger RNA, and LNP Aggregates Using an Ultrawide Pore Size Exclusion Chromatography Column”](#)

Anal. Chem. 2023.




[“Analyzing Encapsulated mRNA with LC, MS, and Calorimetry”](#) Genetic Engineering News, 2021.



[“Liquid Chromatography Methods for Analysis of mRNA Poly\(A\) Tail Length and Heterogeneity”](#)

Anal. Chem, 2023, 95, 38, 14308–14316.



[waters.com/GTx](https://waters.com/GTx)  
[waters.com/GTxResolveSEC](https://waters.com/GTxResolveSEC)  
[waters.com/OligoRPLC](https://waters.com/OligoRPLC)  
[waters.com/OligoStds](https://waters.com/OligoStds)

**Waters™**

Waters, GTXResolve, MaxPeak, XBridge, VanGuard, ACQUITY, BEH, CSH, Gen-Pak, and XSelect are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2024 Waters Corporation. May 24-12739 720008208EN Rev. C

Waters Corporation  
34 Maple Street  
Milford, MA 01757 U.S.A.  
T: 1 508 478 2000  
F: 1 508 872 1990  
[waters.com](https://waters.com)