

# Advances in Downstream PAT for Biologics, Vaccines and Gene Vectors




Daniel Some – Waters | Wyatt Technology  
Contact: [dsome@wyatt.com](mailto:dsome@wyatt.com)



## Abstract

Multi-angle light scattering (MALS) is a versatile technology most commonly found in analytical laboratories for characterizing biophysical properties such as protein or polynucleotide molar mass, size and conformation; size and concentration of virions; and the genomic payload of viral and non-viral gene vectors.

Real-time multi-angle light scattering (RT-MALS) is a novel use of MALS which brings many of those capabilities to preparative systems and process development laboratories. RT-MALS determines key product attributes of biologics, in-line or on-line with downstream and fill-finish processes, and can be used to monitor aggregation, distinguish product from impurities, identify process endpoints and more.



## How MALS works: scattered intensity vs. angle



1. **MW:** the average intensity of scattered light, extrapolated to 0°, with concentration (UV) yields the **weight-average molar mass:**

$$I_{LS}(0) \propto M \times c$$

3. **Particle concentration:** the intensity extrapolated to 0° and particle volume (from the size), yield the **particle concentration N:**


$$I_{particle}(0^\circ) \propto N \times V_{particle}^2$$

2. **Size:** The angular variation of scattered light yields the **z-average rms radius  $R_g$**  (or other dimension).

4. **Payload & titer:** The intensity and 2 UV signals (260 nm & 280 nm) yield capsid and DNA MW, Vg/Cp, empty and full titer.

$$\frac{Vg}{Cp} = \frac{M_{DNA(MALS)}}{M_{DNA(full\ genome)}}$$

## RT-MALS: ultraDAWN™ instrument and OBSERVER™ software



## Case study: Flow-through protein purification

### Flow-through HIC remove mAb aggregates


During the run the column saturates and eventually aggregates break through  
UV does not differentiate monomers, aggregates

#### Goals:

- Identify optimal time to cut pool for maximum yield & specified quality
- Calculate %HMW in final pool

#### RT-MALS solution:

- %HMW (assuming monomer-dimer mixture) is calculated from weight-average molar mass.
- Trigger goes ON when monomers are found (MW > 146 kDa) and OFF when reaching 1.5% dimer.
- RT-MALS optimizes pool regardless of variations in input load and content, buffer variations or column aging, and calculates %HMW in the pool.




Patel et al., mAbs 10(7): 945-950 (2018)

$$\%HMW = \frac{c_d}{c_{tot}} = \frac{M_{w,m-d}}{M_m} - 1$$

$c_d$ : dimer concentration  
 $c_{tot}$ : total protein concentration  
 $M_m$ : monomer molar mass  
 $M_{w,m-d}$ : total solution weight-average molar mass

## Case study: Adenovirus purification



4. **Elution:**
- UV:** indicates "something" is eluting, does not differentiate virus from impurities
  - MALS:** confirms monomeric AdV by constant radius
  - MALS:** Integrates titer peak to determine pooled titer
  - Validation:** Good agreement of size and titer with offline reference data. RT-MALS saves ~ 2 weeks turn-around time for offline analytics

	Size (nm)	Titer (1/mL)
RT-MALS	$R_{geom}$ : 45.7	$7.31 \times 10^{11}$
Reference off-line method	$R_n$ : 48.6	$6.07 \times 10^{11}$

## Case study: Viral vector ultrafiltration


### Ultrafiltration by TFF to increase concentration of drug substance

#### Goals:

- Signal final concentration end point
- Flag possible degradation
- Identify membrane fouling

#### RT-MALS solution:

- OBSERVER monitors trigger condition: titer >  $5 \times 10^8 \text{ mL}^{-1}$
- Response time ~ 12 seconds, does not delay processing
- Monitor average virion size to ensure no degradation
- RT-MALS can also monitor the permeate to identify membrane fouling



## Case study: AAV enrichment

### Ion-exchange chromatography removes empty capsids and other impurities


UV signal does not provide quantitative Vg/Cp or titer, or identify aggregates.

#### Goals:

- Optimize column loading and collect enriched vectors by Vg/Cp in real time;
- Identify aggregates, calculate titer

#### RT-MALS solution:

- Linear gradient (top):** OBSERVER trigger set to collect Vg/Cp > 0.5
- OBSERVER calculates final full and empty titers, Vg/Cp in pool
- Step gradient (bottom):** In Column Load, confirm discard of empty AAVs
- Strip peak identified as aggregates by size  $R_g$



## Conclusions

- RT-MALS - with ultraDAWN and OBSERVER - directly monitors relevant CQAs and other product attributes in downstream processing of biologics like mAbs and viral vectors
- In process development, RT-MALS enhances process understanding and accelerates time to market by cutting down on offline analytics
- In production, RT-MALS increases productivity, yield and quality while adding flexibility for transferring across scales and sites.

