



PFAS Tested Oasis™ SPE Saves Time and Reduces Risk



Oasis WAX/GCB for challenging aqueous samples like surface water, ground water, and wastewater



Oasis GCB/WAX for solid, biosolid, tissue, and food samples



Oasis WAX for drinking water samples and WAX-only clean up methods for other matrices



Quality-control PFAS tested at accredited lab for **reduced risk of contamination**

Reduce time of loose GCB clean-up for complex matrices with dual-phase cartridges

Cartridges meet acceptance criteria for **EPA 1633 and 533**

Description	Part Number
Oasis WAX for PFAS Analysis 6 cc, 500 mg, 30 µm, 30/pk	186009347
Oasis WAX for PFAS Analysis 6 cc, 500 mg, 60 µm, 30/pk	186009346
Oasis WAX for PFAS Analysis 6 cc, 500 mg, 60 µm, 30/pk	186009347
Oasis WAX for PFAS Analysis 6 cc, 150 mg, 30 µm, 30/pk	186009345
Oasis WAX for PFAS Analysis 6 cc, 150 mg, 30 µm, 300/pk	186009344

Description	Part Number
Oasis WAX/GCB for PFAS Analysis 6 cc, 200 mg, 60 µm/50 mg, 30/pk	186011110
Oasis WAX/GCB for PFAS Analysis 6 cc, 200 mg, 60 µm/50 mg, 300/pk	186011111
Oasis GCB/WAX for PFAS Analysis 6 cc, 50 mg/200 mg, 60 µm, 30/pk	186011112
Oasis GCB/WAX for PFAS Analysis 6 cc, 50 mg/200 mg, 60 µm, 300/pk	186011113



Learn More on how you can conquer your PFAS challenges!